

ViscoSert M

Hydroxypropylmethylcellulose 2.0%

Viscoelastic solution for intraocular use



Sterile single-use syringe with 2 ml solution and single-use injection cannula

- ▶ Easy to remove at the end of the surgery
- ▶ large quantity of solution: 2 ml

ViscoSert M

Hydroxypropylmethylcellulose 2.0%

Formula (1 ml)			
Hydroxypropylmethylcellulose	20.00 mg	Sodium acetate	3.90 mg
Sodium chloride	6.40 mg	Sodium citrate	1.70 mg
Potassium chloride	0.75 mg	Water for injection	
Calcium chloride	0.48 mg		
Magnesium chloride	0.30 mg		
Supplied cannula	23 G		

ViscoSert M is distributed by HOYA Surgical Optics GmbH. Manufactured by ALBOMED®, Germany
HOYA Surgical Optics names and logos are registered trademarks of HOYA Surgical Optics, Inc.
©2013 HOYA Surgical Optics, Inc. All rights reserved.

HOYA Surgical Optics GmbH
Headquarter Europe, Middle East and Africa
De-Saint-Exupéry-Str. 8
60549 Frankfurt am Main
Germany

Phone: +49 69 664 268 0
Fax: +49 69 664 268 399
HOYA.com/SurgicalOptics
ifra-info@HOYA.com

Singularly Focused. [Globally Powered.](#)

HOYA
SURGICAL OPTICS