

ViscoSert P

Sodium Hyaluronate 1.0%

Viscoelastic solution for intraocular use



Sterile single-use syringe with 0.85 ml solution and single-use injection cannula

- ▶ High degree of biocompatibility
- ▶ Good space maintenance
- ▶ Controlled capsulorhexis
- ▶ Excellent protection of ocular tissues
- ▶ Easy implantation

ViscoSert P

Sodium Hyaluronate 1.0%

Technical data		Formula (1 ml)	
Sodium Hyaluronate	1.0 %	Sodium hyaluronate	10.0 mg
Molecular Weight (Dalton)	avg. 5 million D	Sodium chloride	8.3 mg
Viscosity	avg. 1 million mPa	Disodium phosphate dihydrate	1.4 mg
Osmolality	310 mOsm/kg	Potassium dihydrogen phosphate	0.26 mg
pH	6.8 – 7.6	Water for injection	q.s.
Storage	2 – 25 °C		
Volume	0.85 ml		
Supplied cannula	27 G		

ViscoSert P is distributed by HOYA Surgical Optics GmbH. Manufactured by Bohus BioTech AB, Sweden
 HOYA Surgical Optics names and logos are registered trademarks of HOYA Surgical Optics, Inc.
 ©2013 HOYA Surgical Optics, Inc. All rights reserved.

HOYA Surgical Optics GmbH
 Headquarter Europe, Middle East and Africa
 De-Saint-Exupéry-Str. 8
 60549 Frankfurt/Main
 Germany

Phone: +49 69 664 268 0
 Fax: +49 69 664 268 399
 HOYA.com/SurgicalOptics
 ifra-info@HOYA.com

Singularly Focused. [Globally Powered.](#)

